



User

Date

Tower

Work report

User

Last name : _____

First name : _____

Phone number : _____

E-mail : _____

Billing address

Laboratory / service : _____

Unit : _____

Supervisor : _____

Address : _____

Country : _____

E-mail : _____

Samples and measurement conditions

N°	Simplified sample identification	Holder	BSML file or specific conditions	Remarks (toxic, PP, etc.)
1				
2				
3				
4				
5				
6				
7				
8				
9				
10				
11				
12				
13				
14				
15				

Sample Holders examples (please ask for more codes) :

SN - Zero Background support

SNF - Floating Zero Background support

PMMA SF (Fine Striations) - SP (Strong Striations) - FP (Flat base)

D8A-MC - Mylar Chamber

Thick Solids (manual mode)

BSML files if applied, using the format :

D8A-T2T-min-max-kV-mA-FF/FV-rot-KEC-Ni-XE-duration

i.e. :

D8A-T2T-01-010-40kV-40mA-FF01-SS-norot-KECM-Ni-SS-XE3-3min

by default, measures are done with 2Theta from 5-70°

Measurement files :

USER-YY-WW-DNN

USER : User Code

Year Week Day Number

example : PROB-19-19-401

for Pierre ROBERT on 09/05/2019